

hyperpigmented skin

Hyperpigmentation can be described as an abnormal excessive pigmentation of the skin, usually on the face. It may be caused by any irritant to the skin, but the main cause is exposure to UV radiation. Nimue's range of products is differentiated as a derma cosmeceutical line for dramatic improvement of hyperpigmented skin.

nimue
SKIN TECHNOLOGY

EXPERIENCE SKIN EXCELLENCE

treatment programme

am

- **Cleansing Gel**
 - Massage skin for 1 minute with damp fingertips
 - Rinse thoroughly
 - Repeat if necessary
 - Dry gently
- **Conditioner**
 - Apply a small amount to the skin, using a damp or dry gauze pad
 - Allow to dry
- **Day Fader**
 - Apply as needed
 - Massage into skin
- **SPF 20**
 - Apply to the skin (over Day Fader) every 2-3 hours
 - Re-apply after exercise or water activities

pm

- **Cleansing Gel**
 - Massage skin for 1 minute with damp fingertips
 - Rinse thoroughly
 - Repeat if necessary
 - Dry gently
- **Conditioner**
 - Apply a small amount to the skin, using a damp or dry gauze pad
 - Allow to dry
- **Exfoliating Enzyme**
 - Apply a thin layer to cover areas treated
 - Leave for 2 minutes
 - Massage with wet fingertips
 - Leave for another 2 minutes
 - Rinse thoroughly
 - Use as prescribed
- **Active Gel/Lotion**
 - Apply a pea-sized drop to the skin
 - A mild stinging sensation may be experienced
 - Do not rinse off
 - Use every second night - gradually increase to daily use
 - Not to be used the night after a clinic treatment
- **Night Fader**
 - Apply as needed
 - Massage into skin

additions



before



niVue skin scanner analysis



after