



# COLD STERILANT®

SAFE COLD STERILISATION & DISINFECTING OF MEDICAL EQUIPMENT

For the safe, cold sterilisation and disinfecting of your surgical instruments and heat sensitive equipment.

## SAFE

Non-toxic, Non-Corrosive, Non-Irritant, Non-Staining, Non-Caustic, Non-Carcinogenic, Biodegradable, pH Neutral. Contains no Aldehydes, Phenol or Chlorine.

## QUICK

Sterilises in 30 minutes and Disinfects in 5 minutes with a residual activity that keeps on killing microbes.

## EASY

Ready to use, no fuss, no smell and no need to rinse.

## ECONOMICAL

Replace only every 30 days after light to moderate use.  
14 days under normal to heavy use. Has a long stable shelf life.

### Cold Sterilant® is available in:

- 1 Litre - Bottle
- 5 Litre - Bulk

